



# HJ Enterprise

## ::: Horizontal High Pressure Cyelndrical Steam Sterilizer :::



### ::: Available Size :::

- 12"x22" Dresing drum 9"x11" - 2 nos.
- 16"x24" Dresing drum 14"x12" - 2 nos.
- 20"x36" Dresing drum 18"x11" - 3 nos.
- 20"x48" Dresing drum 18"x11" - 4 nos.

### ::: Salient Feature :::

- End ring hydraulically formed without joint to sustain Required pressure. All other joints Argon Arc welded.
- Two safety valves.
- Low water level device.
- Pressure lock device.
- Pressure lock device for the door.
- Works on 121°C at 17 psi(1.2kg/cm<sup>2</sup>).
- Hydrostatically tested 3.5kgf/cm<sup>2</sup>.
- Vacuum system incorporated.
- Power: 9KW. 440 V,3-Phase & Neutral.
- First cycle : 90 minutes.
- Subsequent :40 minutes.
- PC /PLC based Auto Control Panel. (optional)
- Vaccuim pump for drying.(optional)

### ::: Other Models :::

- HJE 102 : 18 KW.
- HJE 103 Mini : Regular cum High Speed Best Suited for private nursing homes and all Surgical Branches.
- Boiler/Pipe Steam any combination for all.