



HJ Enterprise

::: Horizontal Rectangular High Pressure Steam Sterilizer :::



::: Salient Features :::

- Best suited for parenteral solutions, injectables, mattresses And hospital dressings.
- Unipoint control station for all phases of the sterilizing cycle In smaller size.
- Operating Temperature 121°C or higher 134°C.
- All necessary safety precautions namely Pressure locking Door, Spring loaded Safety valve, dead weight Safety valve, Etc.. Incorporated.
- Loading carriage and trolley as standard accessories for easy transfer of load.
- Choice of Single Door, Double Door, Recessed or Panel Models.
- Various material combination like Mild Steel, Stainless Steel 304/316 is available.
- Manufacturing as per customers requirement.
- Steam Trap is incorporated in discharge line for fast and perfect elimination of air.

::: Optional Features :::

- Single Door / Double Door.
- Flush Mounted.
- Digital Temperature Recorder.
- Rapid Heating
- Rapid Cooling (Air & Water). High Vacuum Pump.
- PC(Computer)/ PLC based Auto Control Panel.
- Air Oven Steam Pulsing Pressure System

::: Available Sizes :::

600 mm x 600 mm x 1200 mm
600 mm x 900 mm x 1500 mm
900 mm x 900 mm x 1800 mm
1200mm x 1200mm x 2400mm

HJ Enterprise
(Work) : 889/2, G.I.D.C, Makarpura, Baroda (India). 390010.
(Office) :C/1 Gautam Nagar Socitey, Alwa Naka G.I.D.C, Road Manjalpur
Baroda (India) 390011.
Phone : +91-98244 39383, 98257 41041