



HJ Enterprise

::: Vertical Autoclave Double Drum (S.S) H.J.E.-105 :::



::: Specification :::

Double wall unit made of stainless steel chamber and lid made of stainless steel single piece, unit run on power 220V/4000 w 50 Hz double safety valves to ensure the safety of the user and the pressure gauge on the top of the lid, whistle, safety valve also provided on the top of the lid, gasket placed well between the lid and the unit water Lever Indicator to show the level of the water inside the unit, release wave of the steam while the cycle complete, water changing tap, water input funnel and two metal handle to hold the unit while movement, accessories; Heavy duty S.S. Power lead, Silicone rubber gasket, Lead fitting wing key.

HJ Enterprise

(Work) : 889/2, G.I.D.C, Makarpura, Baroda (India). 390010.

(Office) :C/1 Gautam Nagar Socitey, Alwa Naka G.I.D.C, Road Manjalpur
Baroda (India) 390011.

Phone : +91-98244 39383, 98257 41041