



HJ Enterprise

::: Vacuum Type Vertical Sterilizer . H.J.E.- 106 :::



**With Digital Temperature Meter
Steam Trap Is Incorporated In
Discharge Line For Condensed
Water**

- **Steam Working Pressure** 1.26 kgf/cm² (2.2 kgf/cm² in case High speed sterilizer)
- **Steam Working Temperature** 121°C (134°C in case of High speed sterilizer)
- **Hyd. Tested at** Jacket-Twice the working pressure.
Chamber-One & half time the working pressure.
- **Operating Voltage** 400/440, 3 phase, 9000 W, AC supply, 50c/s.
- **Sterilization Period** 20 to 25 minutes(in case of High speed sterilizer 5 to 7 minutes).
- **Steam Exhaust** 5 to 7 minutes(in case of High speed sterilizer 1 minute)

Autoclaves is availabe as per customer's requirement

HJ Enterprise
(Work) : 889/2, G.I.D.C, Makarpura, Baroda (India). 390010.
(Office) :C/1 Gautam Nagar Socitey, Alwa Naka G.I.D.C, Road Manjalpur Baroda (India) 390011.
Phone : +91-98244 39383, 98257 41041